# BT Scan Strip

ZE/BT25  
ZE/BT100

Immunochromatographic test  
for  $\beta$ -lactams and tetracyclines in milk

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## SCOPE AND TEST PRINCIPLE

BT Scan Strip is a qualitative test designed for detection of betalactams ( $\beta$ ) and tetracyclines (T) in milk. This test can be used in raw, skimmed, heated or powder milk from cow, sheep, goat or buffalo. The analysis of raw milk can be performed in individual or commingled milk.

BT Scan Strip is a competitive test with two components: a tube with freeze-dried reagents that has to be mixed with the milk sample, and a strip which is subsequently introduced into the tube. The tube contains receptors and antibodies bound to gold particles and the strips have specific antibodies. If the sample is free from antibiotics  $\beta$  or T, the receptors present in the tube will bind to its specific line of the strip when the liquid runs vertically, and intense red lines will appear. When some antibiotic is present in the sample, it will bind the receptors of the tube, which will not be able to bind to the specific line in the strip, inhibiting the appearance of the red color of the corresponding test line, or decreasing its intensity.

## KIT COMPONENTS

	ZE/BT25	ZE/BT100
INDIVIDUAL TUBES WITH FREEZE-DRIED REAGENT	25	100
STRIP TESTS	25	100
PIPETTE TIPS (200 $\mu$ L)	25	100
PRODUCT CERTIFICATE	1	1
KIT INSTRUCTIONS	1	1

## ADDITIONAL MATERIAL (NOT PROVIDED)

- Micropipette 200  $\mu$ L
- Block incubator at 40°C or FS incubator
- Timer
- Reader IRIS - optional (reference ZE/IRIS)

**SAFETY**

A SAFETY DATA SHEET is available from your local distributor or ZEULAB. This kit should be used following good laboratory practices.

**STABILITY AND STORAGE**

Store the test in dry conditions and under refrigeration (4-12°C) until use. Check the expiry date on the package.

**NOTES**

- The strips are very sensitive to moisture, keep them inside the supplied container and fully closed.
- Let the kit to reach room temperature before its use.
- Milk samples should not be tested if they are older than 2 days.
- A new pipette tip must be used for each sample.
- Read the result within the one minute after incubation. Strips will dry after incubation and colour intensities of lines will sharpen and interpretation will not be valid.
- If the control line does not appear, the assay might have not been run properly or the reagents were damaged. The test must be repeated.
- Make sure that the sample does not reach the area of the test strip arrows. For manipulating the strip, hold it from the opposite end of the arrows area.
- The IRIS reader is available for reading the strips and ensure an objective interpretation of the test results. Results will be available in the cloud.
- Please, contact ZEULAB for further information.

## TEST PROCEDURE

1. Cut the tubes needed, remove the foil and place them in a preheated block incubator at  $40 \pm 1^\circ\text{C}$ . Ensure that the unused tubes remain covered with the foil to prevent humidity.
2. Add 200  $\mu\text{L}$  of milk sample into the tube and mix 3 times using the pipette to dissolve the freeze-dried reagent.
3. Take the strip test and place it into the tube in contact with the mixture of milk and reagent, following the arrows direction.
4. Incubate at  $40 \pm 1^\circ\text{C}$  for 6 minutes.
5. Read the results visually or objectively with IRIS in the following minute after the end of the assay.

## INTERPRETATION OF RESULTS

**Visual reading - See the figure at the end of the document**

NEGATIVE: test line ( $\beta$  or T) is more intense than the control line (C).

POSITIVE: test line ( $\beta$  or T) is as intense as the control line (C) or less visible. The sample tested contains an antibiotic concentration that is equal or higher than the limit of detection. Test line can disappear when sample concentration is well above the limit of detection of the test

INVALID: control line (C) is not shown. This may occur when the assay has not been performed properly, when the matrix interferes or when the reagents are damaged. Repeat the assay with a new strip and contact ZEULAB for further information.

### Objective reading – Reader IRIS

Place the strip in the adapter and couple it with the reader IRIS. Please, check the IRIS reader instructions for further information.

**BT Scan Strip** is an *in vitro* diagnostic kit. In analysis implicating legal processes, the results should be re-evaluated with an official reference method. ZEULAB does not assume any legal responsibility. To the best of our knowledge, the information contained herein is accurate and complete. However, nothing herein shall be construed to imply any warranty or guarantee.

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