Eclipse Farm 4G

ZE/EF4G25

ZE/EF4G50

Tube tests for detection of
inhibitory substances in milk with COMET4

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SCOPE AND PRINCIPLE

ECLIPSE FARM 4G is a qualitative test designed for detection of antibiotics and inhibitors in milk. This test can be used in raw cow milk from individual or commingled milk.

It is based on the inhibition of microbial growth. The tubes provided contain an specific culture medium spread with *Geobacillus stearothermophilus* spores and a pH indicator. When the tubes are incubated in a COMET4 device in absense of antibiotics, spores germinate and cells grow (negative). If samples contain inhibitors at higher concentration than the detection limit, microorganisms will not grow (positive). The system has been optimized by combining the kit with COMET4 equipment, which integrates the processes of incubation, time control, interpretation of results and information in real time.

KIT COMPONENTS

	ZE/EF4G25	ZE/EF4G50
Test tubes	25	50
Mini-pipettes	25	50
Adhesive film	1	2
Kit instructions	✓	✓
Product certificate	✓	✓

ADDITIONAL MATERIAL (NOT PROVIDED)

- COMET4

NOTES

- ✓ Before starting, install and switch on COMET4 device following the manual provided.
- ✓ A new pipette should be used for each sample. To use the provided pipettes: squeeze the bulb, put the stem in the sample and release the bulb to draw up liquid. Excess liquid is siphoned automatically into overflow tanks. Squeeze the bulb again to release 100 µl.
- ✓ This test is extremely sensitive to antibiotics and other antibacterial substances such as detergents and disinfectants. Any contamination with those substances should be prevented.
- ✓ Milk contains natural inhibitors which do not interfere with the result. However, these substances appear in high concentrations in colostrum, at the end of lactation and in mastitis cases, and they may alter the test results.

TEST PROCEDURE (Flowchart Procedure on the last page)

Read the Notes section carefully.

1. Switch COMET4 on and follow the instructions provided with the device. Wait until the optimum temperature is reached.
2. Cut an Eclipse Farm 4G tube with scissors (take care not to damage the foil on the remaining tubes).
3. Add 100 µL of milk with the provided pipette and seal the tube with the adhesive film.
4. Place the tube into the chosen position of COMET4 and press "Run Test" on the App.
5. Result will be notified in real time. Check the results as indicated in COMET4 instructions.

STABILITY AND STORAGE

The kit components should be stored at 4-12°C and protected from light. The kit is provided with a minimum shelf life of 9 months. See the kit expiry date on the package.

SAFETY

Good laboratory practices are recommended when using this kit. A Safety Data Sheet (SDS) is available from your local distributor or ZEULAB upon request.

ECLIPSE is an *in vitro* diagnostic kit for antibiotics screening in milk. Analyses implicating legal processes should be carried out in duplicate or triplicate and confirmed by an official reference method. ZEULAB does not assume any legal responsibility.

To the best of our knowledge, the information contained herein is accurate and complete. However, nothing herein shall be construed to imply any warranty or guarantee.

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